

You want results,
it's that simple.
Nothing delivers
results like
Chromolite™.

Safe and Effective
for Men and Women

PERMANENT HAIR REDUCTION*

Treats the removal of unwanted hair in the following areas of the body:

- Face
- Neck
- Back
- Legs
- Bikini Area
- Under Arm
- Arms
- Hands
- Feet

PHOTO REJUVENATION

Treats Benign Pigmented Epidermal Lesions Including:

- Dyschromia
 - Hyperpigmentation
 - Melasma
 - Ephelides (freckles)
- Cutaneous Lesions Including:
 - Warts
 - Scars
 - Striae

Benign Cutaneous Vascular Lesions Including:

- Port Wine Stains
- Hemangiomas
- Facial, Truncal and Leg Telangiectasias
- Erythema of Rosacea
- Angiomas
- Spider Angiomas
- Poikiloderma of Civatte
- Leg Veins
- Venous Malformations

ACTIVE ACNE

Treats:

- Moderate Inflammatory Acne Vulgaris

* Permanent Hair Reduction is defined as a long-term stable reduction in the number of hairs re-growing after a treatment regimen.

Brought to you by

Genesis
Biosystems

405 State Hwy. 121 Bypass | Suite #D-120 | Lewisville, Texas 75067
Toll Free: (888) 577-7335 | Phone: (972) 315-7888
Fax: (972) 315-7818 | www.chromolite.us

MA1-0179-004A

chromolite



Achieve
A Better
Appearance
At The Speed
of Light!

INTENSE PULSED LIGHT TREATMENT FOR ACNE,
PERMANENT HAIR REDUCTION & SKIN REJUVENATION

Advanced technology.

Quicker treatment. **Better results.**



You want results. It's that simple. Nothing delivers results like Chromolite™.

Chromolite™ is a state-of-the-art hair reduction and skin rejuvenation system specially designed for the treatment of unwanted hair, photo rejuvenation, and acne for use with skin types I-V.

Chromolite™ uses proprietary Smartlite™ treatment technology which produces high levels of light energy –and fast results.

Our Smartlite™ Flash Lamp hand piece efficiently covers large treatment area with little light energy waste.

The result is the highest quality treatment in the quickest time.

Generally, a course of four to six treatments is required to achieve significant results.



PERMANENT HAIR REDUCTION

Chromolite's™ Smartlite™ delivers intense pulsed light to increase the efficiency of results for hair reduction. Research has shown that in hair reduction, it is important to destroy the germ cells (bulge), as well as the hair root (bulb). Smartlite's™ yellow light targets the blood vessels around the bulge and its red light spectrum targets the melanin within the hair shaft. This leads to destruction of both the bulge and the bulb for the most effective possible treatment.



ACTIVE ACNE

Chromolite™ uses the yellow section of its light to improve acne because this light level reacts with porphyrin to produce oxygen which destroys the unwanted bacteria in the patient's skin.

Chromolite™ results are improved over other light-based systems

Chromolite™ results are improved over other light-based systems, because Chromolite's™ yellow light acts to produce an anti-inflammatory response in the skin as well.



Chromolite™ uses Smartlite™ technology to achieve enhanced results

PHOTO REJUVENATION

Photo rejuvenation is intense pulsed light applied to the skin to improve the overall appearance. In lighter skin types, Chromolite's™ yellow light targets the blood vessels in the skin causing them to shrink while its red light targets any unwanted pigmentation and stimulates collagen growth. In darker skin types, the yellow light spectrum is filtered out. The result is a smoother complexion with even skin tone.



chromolite 

www.chromolite.us

